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Customized PEEK Implant

Intended Use

Ad Mirabiles Customized PEEK implants are used to treat osteotomies or augmentations of bony defects. The implants may be used in combination with autograft or allograft bone tissue but not with other medical devices except those required to fixate the implant. The implants may be used in trauma medicine, as a cranial fastener in neurosurgery procedures, in plastic surgery or after resection of a bone tumour. Ad Mirabiles Customized PEEK implants are applicable in cases when standard solutions are not available or offer less therapeutic value for the patient. The implant is intended exclusively for the European market (including EFTA states) and is not to be used in the US or Canada.

Users

These products must be used only by medical professionals with the qualifications necessary for these procedures who are familiar with the specific clinical case for which the intended implant has been produced.

Intended use of the product as defined by the manufacturer:

Patient-specific cranial and craniomaxillofacial Ad Mirabiles implants are intended to fill bony defects or correct the patient's anatomy. The design of an Ad Mirabiles Customized PEEK implant is based on the patient's individualised CT scan and thus adapted to the specific anatomy of the patient.

Contraindication

- Existing or suspected infections close to the site of implant.
- Insufficient bone density or quality for implant attachment.
- Patients with poor cooperation during rehabilitation, for example due to drug abuse, alcoholism or for mental reasons.
- Treatment with an implant is not recommended if the patient belongs to a risk group such as haemophilia, unsatisfactory wound healing or reduced immunity.

Influencing factors

Occupational and/or physical activities that are associated with substantial forces and hence stress on the body are risky. In extreme cases, these activities may result in product failure and thus may even cancel out the actual purpose of the intervention.

In rare cases, a degenerative disorder at the time of the implantation may be so far advanced that the expected lifetime of the implant may be significantly reduced. In such cases, the products should simply be considered as measures to delay the condition or provide temporary relief.

If there is justified cause to suspect a sensitivity or allergy to the materials used in the implant, appropriate tests should be undertaken by the attending physician prior to selecting materials and performing the implantation.

Customized PEEK Implant

Adverse effects

In most cases, secondary reactions not directly related to the implant design itself or to the material may cause ongoing complications. These may include:

- Loosening of the implant due to inadequate attachment
- Bone necrosis, osteoporosis, insufficient revascularisation, bone resorption and poor bone formation that may cause premature loss of fixation.
- Irritation of the soft tissue and/or neural damage through surgical trauma.
- Infections occurring early or late, both on the surface or within the tissue.
- Increased fibrotic tissue reactions in the region surrounding the surgical intervention site.
- Complications at the time of removal of the implant (if required) due to improper explanation.

Preparation before implanting the Ad Mirabiles Customized PEEK Implant

Inspection

Implants should be inspected for damage such as scratches, dents, wear, discolouration or residues. Damaged implants must be returned immediately to the manufacturer. Implants should not come into contact with foreign substances other than those used for disinfection, cleaning or sterilisation (see the section below).

Customized PEEK Implant

Preliminary cleaning

Ad Mirabiles Customized PEEK implants are delivered pre-cleaned. Prior to sterilisation, Ad Mirabiles implants must undergo a preliminary cleaning according to the guidelines of EN ISO 17664.

Cleaning: Automated/mechanical washing method

Equipment: Washer machine/disinfector, enzymatic cleaner or detergent solution.

Use the following cycle parameters:

Cycle	Minimum time (minutes)	Minimum water temperature	Cleaning agent type
Pre-wash	2	Cold tap water	N/A
Wash I	2	Cold tap water (<40°C)	Cleaning agent* pH 7-9.5
Wash II	5	Warm tap water (>40°C)	Cleaning agent * pH 7-9.5
Rinse	2	Warm, deionised or demineralised water (> 40°C)	N/A
Thermal disinfection	5	> 93°C	N/A
Drying	40	> 90°C	N/A

*Please consult the documentation provided by the manufacturer

Packaging

Prior to sterilisation, cleaned and dried implants should be placed in appropriate sterilisation sleeves or in a reusable solid sterilisation vessel, for example a sterile barrier system according to ISO 11607. The implants and sterile barrier system must be protected from contact with pointed or sharp objects which may damage the surface of the product or the sterile barrier system.

Customized PEEK Implant

Sterilisation

All Ad Mirabiles Customized PEEK implants are delivered in a **NON-STERILE** state! Prior to use, the implants must be sterilised according to the following parameters in accordance with EN ISO-17664:

Sterilisation method (cycle)	Sterilisation time	Sterilisation temperature	Drying time
Steam sterilisation (fractionated prevacuum) (minimum 3 intervals)	min. 4 minutes	min. 132°C max. 138°C	20-60 minutes
Steam sterilisation (fractionated prevacuum) (minimum 3 intervals)	min. 3 minutes	min. 134°C max. 138°C	20-60 minutes

Please follow the operating instructions of the autoclave manufacturer and the recommended guidelines for maximum sterilisation pressure. The autoclave must be properly installed, maintained and calibrated. End users should use only approved sterilisation barriers (for example wraps, bags or containers) to package sterilised end products.

The user bears ultimate responsibility for validation of sterilisation techniques and equipment.

Products must be stored in a dry and dust-free environment after sterilisation.

Distributors and manufacturer assume no responsibility for customer sterilisation methods that do not comply with the recommendations and standards named above.

Additional information

For additional information, please consult the respective national laws and guidelines. Please follow the in-house guidelines and procedural instructions of the clinic and the recommendations and instructions of the manufacturer of the cleaning and disinfectant agents and the reprocessing system.

Customized PEEK Implant

Important

The instructions described here for clinical preparation have been reviewed by Ad Mirabiles and are suitable for preparing non-sterile medical products. Nevertheless, the user remains responsible for ensuring that the preparation is undertaken by suitably trained personnel using appropriate and correctly installed, maintained and tested/validated systems and materials in order to achieve the desired result. Any deviation from the instructions described above should be properly checked for efficacy and potential adverse effects.

Please note

- Ad Mirabiles Customized PEEK implants must be used only by users who have the qualifications necessary for this type of intervention and who are familiar with the clinical case.
- As the manufacturer, Ad Mirabiles recommends that the user read all available documentation prior to first use and contact other users who have experience with this type of treatment.
- All Customized PEEK implants made by Ad Mirabiles are intended for single use in one patient and must never be reused (**SINGLE USE!**).
- Ad Mirabiles Customized PEEK implants have been especially developed and produced for the specific anatomy of the patient.
- Customized PEEK Implants made by Ad Mirabiles using PEEK are attached using standard cranial and craniofacial plates to the native bone in the patient. The holes for the osteosynthesis products (plates, screws etc.) must be drilled into the implant by the surgeon prior to implantation.
- For attaching the PEEK implants to a patient's cranial bone, approved osteosynthesis systems (plates and screws) made by various manufacturers are permitted and must be used according to their instructions for use.
- Prior to using the plate-screw system, the user should clarify with the manufacturer whether the chosen system is approved for this indication.
- Customized PEEK implants are delivered in a non-sterile state and have 2.0 mm drainage holes.

The following preparation steps must be carried out outside of the surgical area

- Adaptation (e.g. drilling, milling etc.) of the implants.
- Drillings holes in the PEEK implant for preparation of the osteosynthesis materials to be used to attach the implant to bone.
- Attaching osteosynthesis material to the implant to attach it to bone (see specific surgical technique for plates and screws from the manufacturer).
- Implants should be kept covered until they are ready to use to avoid contamination.

Single use products

Products intended for single use only must not be reused (see "Explanation of symbols").

Customized PEEK Implant

Ad Mirabiles advises against clinical reprocessing of contaminated implants. Products contaminated with blood, tissue and/or bodily fluids and substances should never be reused and must be disposed of according to the guidelines and regulations of the hospital. Even if components appear to be intact after use, minor defects and invisible material faults may cause material fatigue and hence failure of the implant.

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Do not use product if package is damaged